



# MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI

BOLOGNA, ROYAL HOTEL CARLTON  
12 OTTOBRE 2026

## Pitfalls nella diagnostica citofluorimetrica

*Prof. Francesco Buccisano*

## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI

### Disclosure: Francesco Buccisano

| Company name           | Research support | Employee | Consultant | Stockholder | Speakers bureau | Advisory board | Other |
|------------------------|------------------|----------|------------|-------------|-----------------|----------------|-------|
| JAZZ<br>PHAMACEUTICALS |                  |          |            |             | X               | X              |       |
| NOVARTIS               |                  |          |            |             | X               |                |       |
| MENARINI<br>STEMLINE   |                  |          |            |             | X               |                |       |
| ASTELLAS               |                  |          |            |             | X               |                |       |
| SERVIER                |                  |          |            |             | X               |                |       |

BOLOGNA, 12 OTTOBRE 2026

## Pitfalls in MFC diagnostics of BPDCM and SM

Historical rarity of disease, Many nomenclature changes (WHO 2008, WHO 2016)

Difficulty in clinico-pathologic recognition/diagnosis

Morphology not characteristic, unspecified reason for testing

Minority population, fibrosis, hemodilution

Concomitant hematological disease (~20%)

Unmet technical requirement and gap of knowledge

## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI

### Diagnostic/clinical «awareness» of BPDCN

|                                                                                                                           | Replies                  | Percentages |
|---------------------------------------------------------------------------------------------------------------------------|--------------------------|-------------|
| Do you have a BPDCN immunophenotypic panel (CD4-CD56-CD123-BDCA2/CD303) that you use on a regular basis?                  | ✓ Yes                    | 33          |
|                                                                                                                           | ✓ No                     | 17          |
|                                                                                                                           | ✓ Only in selected cases | 50          |
| Your immunophenotypic panel includes markers that exclude BPDCN, such as: CD3, CD19, myeloperoxidase, CD11c, CD14, CD64?  | ✓ Yes                    | 0           |
|                                                                                                                           | ✓ No                     | 100         |
| Does your histochemical panel include CD4, CD56, CD123, BDCA2, TCL1, or excludes: MPO, CD3, CD20 o CD19, CD163, Lysozima? | ✓ Yes                    | 57          |
|                                                                                                                           | ✓ No                     | 43          |
| Do you routinely perform CD123 in acute myeloid leukemia at diagnosis?                                                    | ✓ Yes                    | 59          |
|                                                                                                                           | ✓ No                     | 41          |
| Do you routinely perform a cerebral fluid exam at diagnosis?                                                              | ✓ Yes                    | 70          |
|                                                                                                                           | ✓ No                     | 30          |
| Do you use a different therapeutic approach between BPDCN patients with leukemic spread or without?                       | ✓ Yes                    | 26          |
|                                                                                                                           | ✓ No                     | 74          |
| What type of chemotherapy do you use? AML-oriented, ALL-oriented or NHL-oriented?                                         | ✓ AML-like               | 39          |
|                                                                                                                           | ✓ ALL-like               | 44          |
|                                                                                                                           | ✓ NHL-like               | 17          |


NHL, non-Hodgkin lymphoma.

***41% do not perform routinely CD123 at diagnosis, 50% test BPDCN only in selected cases***

Valentini C et al., Blood Adv 2021 (letter)

BOLOGNA, 12 OTTOBRE 2026

## The role of MFC in BPDCN



## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI

### BPDCN diagnosis: multiparameter flow-cytometry

| ANTIGENS THAT HELP CONFIRM THE DIAGNOSIS OF BPDCN <sup>1-4**§</sup>                |                                                                        |
|------------------------------------------------------------------------------------|------------------------------------------------------------------------|
| CD123+<br>CD4+<br>CD56+<br>TCL1+<br>TCF4+<br>CD303+ (BDCA2+)                       | CD3-<br>CD19-<br>TIA1-<br>Granzyme B-<br>Myeloperoxidase-              |
| ANTIGENS THAT EXCLUDE THE DIAGNOSIS OF BPDCN <sup>1</sup>                          |                                                                        |
| <b>Pan T-cell markers</b><br>CD3<br><br><b>Pan B-cell markers</b><br>CD20<br>CD79a | <b>Myeloid markers</b><br>CD11c<br>CD34<br>Lysozyme<br>Myeloperoxidase |

CD123, CD4, and CD56, completes a signature marker triad that can differentiate BPDCN from other hematologic malignancies<sup>4,5§</sup>

BPDCN is characterized by the expression of the pDC-associated marker, CD123<sup>3§</sup>

#### Immunophenotypic diagnostic criteria:

Expression of CD123 and one other pDC marker (TCF4, TCL1, CD303) in addition to CD4 and/or CD56, or,  
Expression of any three pDC markers and absent expression of all expected negative markers

## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI

# MFC panel construction (Euroflow approach)




**Table 14.** The EuroFlow AML/MDS antibody panel<sup>a,b</sup>

| Tube | PacB  | PacO | FITC           | PE     | PerCPy5.5 | PECy7 | APC            | APCH7 | Aim                                                                                                      |
|------|-------|------|----------------|--------|-----------|-------|----------------|-------|----------------------------------------------------------------------------------------------------------|
| 1    | HLADR | CD45 | CD16           | CD13   | CD34      | CD117 | CD11b          | CD10  | Diagnosis and classification, neutrophilic maturation, PNH                                               |
| 2    | HLADR | CD45 | CD35           | CD64   | CD34      | CD117 | CD300e (IREM2) | CD14  | Diagnosis and classification, monocytic maturation, PNH                                                  |
| 3    | HLADR | CD45 | CD36           | CD105  | CD34      | CD117 | CD33           | CD71  | Diagnosis and classification, erythroid maturation                                                       |
| 4    | HLADR | CD45 | NuTdT          | CD56   | CD34      | CD117 | CD7            | CD19  | Aberrant expression of lymphoid markers, abnormal B lymphoid maturation                                  |
| 5    | HLADR | CD45 | CD15           | NG2    | CD34      | CD117 | CD22           | CD38  | Aberrant marker expression, stem cells                                                                   |
| 6    | HLADR | CD45 | CD42a and CD61 | CD203c | CD34      | CD117 | CD123          | CD4   | Diagnosis and classification of AML Megakaryocytic, basophilic, and plasmacytoid dendritic cell lineages |
| 7    | HLADR | CD45 | CD41           | CD25   | CD34      | CD117 | CD42b          | CD9   | Characterization of megakaryoblastic leukemia, and systemic mastocytosis                                 |

Abbreviations: AML, acute myeloid leukemia; APC, allophycocyanin; BB, backbone; BM, bone marrow; Cy7, cyanin7; FITC, fluorescein isothiocyanate; H7, hiline7; MDS, myelodysplastic syndrome; Nu, nuclear; PacB, pacific blue; PacO, pacific orange; PE, phycoerythrin; PerCPy5.5, peridinin-chlorophyll-protein-cyanin5.5; PNH, paroxysmal nocturnal hemoglobinuria. <sup>a</sup>Further information about markers and hybridomas is provided in the Appendix. <sup>b</sup>A total of 96 BM samples were evaluated for selection of the BB markers. An additional 84 BM AML samples were evaluated with this version (final) of the panel.





| AML vs BCP-ALL |     | AML vs T-ALL |     | BCP-ALL vs T-ALL |     |
|----------------|-----|--------------|-----|------------------|-----|
| CD19           | 32% | CyMPO        | 26% | CD19             | 23% |
| CyMPO          | 27% | CD7          | 26% | CyCD3            | 22% |
| CyCD79a        | 20% | CyCD3        | 26% | CD7              | 22% |
| CD45           | 8%  | CD34         | 15% | CyCD79a          | 14% |
| CyCD3          | 6%  | SmCD3        | 6%  | CD34             | 9%  |

|   |       |      |                |        |      |       |       |     |                                                                                                          |
|---|-------|------|----------------|--------|------|-------|-------|-----|----------------------------------------------------------------------------------------------------------|
| 6 | HLADR | CD45 | CD42a and CD61 | CD203c | CD34 | CD117 | CD123 | CD4 | Diagnosis and classification of AML Megakaryocytic, basophilic, and plasmacytoid dendritic cell lineages |
|---|-------|------|----------------|--------|------|-------|-------|-----|----------------------------------------------------------------------------------------------------------|

## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI

### MFC panel construction (Euroflow approach)



**CD45/CD42A-CD61/CD203C/CD34/CD117/CD123/CD4**


Van Dongen JJ et al., Leukemia 2012; Martin-Martín L et al., Oncotarget 2015

BOLOGNA, 12 OTTOBRE 2026



## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI

### MFC panel construction (MD Anderson proposal)





|          | Antibody list                                              |
|----------|------------------------------------------------------------|
| Panel #1 | CD7/CD33/CD19/CD34/CD13/CD2/CD38/CD45                      |
|          | HLA-DR/CD117/CD4/CD34/CD123/CD38/CD45                      |
|          | CD41/CD36/CD56/CD34/CD64/HLA-DR/CD14/CD45                  |
|          | CD5/CD25/CD22/CD34/CD38/CD15/CD45                          |
|          | TdT/MPO/CD34/CD3/cytoCD3/CD45 (cyto tube)                  |
| Panel #2 | HLA-DR/CD64/CD4/CD33/CD56/CD45<br>CD303/CD123              |
| Panel #3 | CD2+/CD4/CD7-/CD38+/CD45/CD56/<br>CD64/CD123/CD303/ HLA-DR |



## Analytical approach to BPDCN






## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI






**BPDCN**

BOLOGNA, 12 OTTOBRE 2026

## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI






Bright  
CD64  
positivity









Bright  
CD33  
positivity

**pCD-like AML**

## Different degrees of CD123 expression






Caso clinico basato sull'esperienza personale del medico

CD34 / CD304 / HLA-DR / CD64 / CD4 / CD33 / CD56 / CD45 / CD303 / CD123

BOLOGNA, 12 OTTOBRE 2026


# BPDCN Is Often Mistaken for Other Hematologic Malignancies




1. Goyal A, et al. In: Carter JB, et al. eds. Atlas of Cutaneous Lymphomas: Classification and Differential Diagnosis. Cham, Switzerland: Springer International; 2015 2. Bueno C, et al. Haematologica. 2004 3. Pagano L, et al. Br J Haematol. 2016

### pDC-AML: specific characteristics

- ✓ AML with  $\geq 2\%$  plasmacytoid dendritic cells (pDC) has been recently described as pDC-AML
- ✓ Low incidence ( $\sim 5\%$  of AML)
- ✓ pDC expansion with frequent RUNX1 mutations
- ✓ Poor clinical outcome



## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI




1. Immature BPDCN (Cyan):  
 $CD34^{POS} CD123^{DIM}$
2. Intermediate BPDCN (light blue):  $CD34^{BRIGHT} CD123^{BRIGHT}$
3. Mature BPDCN (blue):  
 $CD34^{NEG} CD123^{BRIGHT}$

## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI

### Sensitivity of MFC in detection of abnormal mastcells

- 1 million cells or all sample in the tube at a median flow rate.
- Mast cells were gated based on CD117(bright+)/side scatter (A), followed by CD117+CD45+ gating (B)
- Nonspecific binding and no overlap with myeloblasts were further ensured by examining the CD117/HLA-DR plot (C).
- The mast cells were then analyzed for the expression of CD2, CD25, and CD30 (D).
- An expression of 10% or more positive cells compared with control was considered positive for that marker



**CD2, CD25, CD30, CD45,  
CD117, HLA-DR**

|                      | All SM, No. (%) (n = 87) <sup>a</sup> |                    |                | Reactive, No. (%) (n = 40) |                    |                | P value |
|----------------------|---------------------------------------|--------------------|----------------|----------------------------|--------------------|----------------|---------|
|                      | Uniformly positive                    | Partially positive | Positive total | Uniformly positive         | Partially positive | Positive total |         |
| CD2                  | 41 (47)                               | 27 (31)            | 68 (78)        | 0 (0)                      | 0 (0%)             | 0 (0)          | <.0001  |
| CD25                 | 83 (95)                               | 2 (2)              | 85 (98)        | 3 (8)                      | 2 (5)              | 5 (13)         | <.0001  |
| CD30                 | 48 (55)                               | 30 (34)            | 78 (90)        | 0 (0)                      | 5 (13)             | 5 (13)         | <.0001  |
| 1 abnormality        | —                                     | —                  | 87 (100)       | —                          | —                  | 10 (25)        | <.0001  |
| 2 or 3 abnormalities | —                                     | —                  | 80 (92)        | —                          | —                  | 0 (0)          | <.0001  |

# Improved diagnostic screening and classification of clonal mast cell diseases by ultrasensitive KIT p.D816V detection

These results confirm the high specificity and ultra sensitivity of the Flow-SuperRCA assay for the detection of *KIT* p.D816V, emerging as a well-suited whole-blood and whole-BM test for diagnostic screening and classification of patients with MCAS and mastocytosis


### Context of Research

- More than 90% of patients with **systemic mastocytosis** have *KIT* p.D816V-positive cells in their bone marrow (BM). However, this mutation often remains unidentified in blood samples from cases with a low mutant cell burden.

### Methods


- **Patients:** A total of 337 adult patients with **mast cell activation syndromes (MCAS)** (n = 117) and mastocytosis (n = 220) plus 98 controls

- ***KIT* p.D816V detection by FlowSuperRCA** (vs ASOqPCR)




### Findings

Enhanced *KIT* p.D816V detection through Flow-SuperRCA







# Clinical presentation







- Referral to an experienced center with expert pathology and clinical trials (preferred)  
▼
- Evaluation:
  - Dermatology with skin biopsy
  - Hematology
    - Bone marrow aspirate/biopsy
    - Flow-cytometry, including CD123, CD4, CD56
    - Karyotype/cytogenetics
    - Myeloid panel DNA-seq
  - PET/CT scan
  - Lumbar puncture with prophylactic IT chemo

## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI







BM







CSF

## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI

DIAGNOSIS







MRD



BOLOGNA, 12 OTTOBRE 2026

## MASTOCITOSI E BPDCN: LE DIAGNOSI DIFFICILI

**BPDCN**  
**cutaneous**  
**infiltrate**



## CONCLUSIONS

- ✓ BPDCN and SM are rare diseases that require a multidisciplinary team to allow a timely diagnosis and staging of disease
- ✓ The hematology laboratory has a key role in two ways:
  - Suggest a possible diagnosis in the context of a leukemia screening procedure
  - Confirm diagnosis and organ involvement during the staging procedure
- ✓ Unmet reagent setting and gap of knowledge need to be addressed